

Anti-Nadph Oxidase 4 (Nox4) Polyclonal Antibody For Wb, Ihc-P, Elisa - Q9Nph5

Item Number: CM0018290

Introduction

KO Validated rabbit polyclonal antibody against NADPH oxidase 4 (NOX4), suitable for Western blot, immunohistochemistry (paraffin), and ELISA. Cross-reacts with human, mouse, and rat NOX4. Ideal for redox signaling and oxidative stress studies.

[Learn More](#)

Parameter	Value
Product Name	[KO Validated] NADPH oxidase 4 (NOX4) Rabbit pAb
Remarks/Alias	KOX; KOX-1; RENOX; 4)
Species	Human
GeneID (Human)	50507
GeneID	50507
Immunogen	Synthetic peptide A synthetic peptide corresponding to a sequence within amino acids 328-578 of human NADPH oxidase 4 (NOX4) (NP_058627.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9NPH5
Protein Weight	67kDa
Shipping	Ice bag