

Anti-Nadph Oxidase 4 (Nox4) Rabbit Polyclonal Antibody For Wb, Ihc-P, Elisa - Q9Nph5

Item Number: CM0021420

Introduction

Rabbit polyclonal antibody targeting NADPH oxidase 4 (NOX4), validated for WB, IHC-P, ELISA in human, mouse, rat. For studying oxidative stress, redox signaling, and kidney function. Recombinant immunogen (aa 328-578).

[Learn More](#)

Parameter	Value
Product Name	NADPH oxidase 4 (NOX4) Rabbit pAb
Remarks/Alias	KOX; KOX-1; RENOX; NADPH oxidase 4 (NOX4)
Species	Human
GeneID (Human)	50507
GeneID	50507
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 328-578 of human NADPH oxidase 4 (NOX4) (NP_058627.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9NPH5
Protein Weight	67kDa
Shipping	Ice bag