


Anti-Nae1 Rabbit Polyclonal Antibody For Wb And Elisa - Q13564

Item Number: CM0030585



Introduction

High-quality NAE1 rabbit polyclonal antibody validated for WB and ELISA. Reacts with human, mouse, rat. Target is NEDD8-activating enzyme E1 regulatory subunit, involved in neddylation and cell cycle.

[Learn More](#)

Parameter	Value
Product Name	NAE1 Rabbit pAb
Remarks/Alias	HPP1; ula-1; APPBP1; NEDFIH; A-116A10.1; NAE1
Species	Human
GeneID	8883 (Human)
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-220 of human NAE1 (NP_003896.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q13564
Protein Weight	60kDa
Shipping	Ice bag