

Anti-Naga Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - P17050

Item Number: CM0026694



Introduction

High-quality anti-NAGA rabbit polyclonal antibody validated for WB, IF/ICC, and ELISA. Targets alpha-N-acetylgalactosaminidase, a lysosomal enzyme critical for glycolipid metabolism. Cross-reacts with human, mouse, and rat.

[Learn More](#)

Parameter	Value
Product Name	NAGA Rabbit pAb
Remarks/Alias	GALB; D225674; NAGA
Species	Human
GeneID (human)	4668
GeneID	4668
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 312-411 of human NAGA (NP_000253.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P17050
Protein Weight	47kDa
Shipping	Ice bag