

Anti-Nags Rabbit Polyclonal Antibody For Wb, Elisa - Q8N159

Item Number: CM0014392



Introduction

NAGS rabbit polyclonal antibody validated for WB and ELISA. Recognizes human, mouse, and rat N-acetylglutamate synthase, a mitochondrial urea cycle enzyme. Ideal for metabolic disorder research. - Q8N159

[Learn More](#)

Parameter	Value
Product Name	NAGS Rabbit pAb
Remarks/Alias	AGAS; ARGA; NAGS
Species	Human
GeneID (human)	162417
GeneID	162417
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 295-534 of human NAGS (NP_694551.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q8N159
Protein Weight	58kDa
Shipping	Ice bag