

Anti-Nags Rabbit Polyclonal Antibody - Q8N159

Item Number: CM0020472



Introduction

High-quality rabbit polyclonal antibody against human N-acetylglutamate synthase (NAGS), validated for WB and ELISA. Cross-reacts with human and mouse samples. Ideal for urea cycle research and diagnostic assay development.

[Learn More](#)

Parameter	Value
Product Name	NAGS Rabbit pAb
Remarks/Alias	NAGS; AGAS; ARGAs; N-acetylglutamate synthase
Species	Human
Gene ID	162417
Immunogen	Recombinant protein of human NAGS
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q8N159
Protein Weight	58kDa
Shipping	Ice bag