


Anti-Naip Rabbit Monoclonal Antibody For Ihc-P, Elisa - Q13075

Item Number: CM0002357



Introduction

High-quality anti-NAIP rabbit monoclonal antibody validated for IHC-P and ELISA. Reacts with human and mouse NAIP (BIRC1). Ideal for apoptosis and inflammasome research.

[Learn More](#)

Parameter	Value
Product Name	NAIP Rabbit mAb
Remarks/Alias	BIRC1; NLRB1; psiNAIP; NAIP
Species	Human
GeneID (Human)	4671
GeneID	4671
Immunogen	Synthetic Peptide A synthetic peptide corresponding to a sequence within amino acids 1-100 of human NAIP (Q13075)
Source	Rabbit
Category	Monoclonal Antibodies
Application	IHC-P, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q13075
Protein Weight	160kDa
Shipping	Ice bag