

Anti-Naip Rabbit Polyclonal Antibody For Wb, Elisa - Q13075

Item Number: CM0020912



Introduction

Rabbit polyclonal antibody against human NAIP (BIRC1/NLRB1) validated for WB and ELISA. Detects apoptosis inhibitor and innate immune sensor, with rat cross-reactivity. Ideal for neuronal apoptosis and inflammasome research.

[Learn More](#)

Parameter	Value
Product Name	NAIP Rabbit pAb
Remarks/Alias	BIRC1; NLRB1; psiNAIP; NAIP
Species	Human
GeneID (Human)	4671
GeneID	4671
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 750-850 of human NAIP (NP_075043.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Rat
SWISS	Q13075
Protein Weight	160kDa
Shipping	Ice bag