


Anti-Na⁺/K⁺-ATPase Rabbit Polyclonal Antibody For Wb And Elisa - P05023 / P50993 / P13637 / Q13733

Item Number: CM0022122



Introduction

Rabbit polyclonal antibody targeting human Na⁺/K⁺-ATPase alpha-1 subunit (ATP1A1). Validated for WB and ELISA, cross-reacts with rat. Ideal for ion transport and membrane potential studies.

[Learn More](#)

Parameter	Value
Product Name	Na ⁺ /K ⁺ -ATPase Rabbit pAb
Remarks/Alias	ATP1A1, ATPase Na ⁺ /K ⁺ transporting subunit alpha 1, Na ⁺ /K ⁺ -ATPase
Species	Human
GeneID	476/477/478/480
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 386-595 of human Na ⁺ /K ⁺ -ATPase (NP_000692.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Rat
SWISS	P05023/P50993/P13637/Q13733
Protein Weight	74kDa/109kDa/112kDa/113kDa
Shipping	Ice bag