


Anti-Na⁺/K⁺-ATPase Rabbit Polyclonal Antibody For Wb - P05023

Item Number: CM0026181



Introduction

Na⁺/K⁺-ATPase Rabbit Polyclonal Antibody, validated for WB, IF-P, IHC-P, ELISA in human, mouse, rat. Recombinant immunogen spanning aa 1-100 of human ATP1A1. Ideal for studying ion transport, membrane potential, and neurological disorders.

[Learn More](#)

Parameter	Value
Product Name	Na ⁺ /K ⁺ -ATPase Rabbit pAb
Remarks/Alias	CMT2DD; HOMGSMR2; Na ⁺ /K ⁺ -ATPase
Species	Human
GeneID (Human)	476
GeneID	476
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-100 of human Na ⁺ /K ⁺ -ATPase (NP_000692.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF-P, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P05023
Protein Weight	113kDa
Shipping	Ice bag