

Anti-Nanog Monoclonal Antibody For Wb, If-P, Elisa - Q9H9S0

Item Number: CM0006945



Introduction

Rabbit monoclonal antibody against human Nanog, validated for WB, IF-P, ELISA. Reacts with human, mouse, rat. Suitable for stem cell and developmental biology research.

[Learn More](#)

Parameter	Value
Product Name	Nanog Rabbit mAb
Remarks/Alias	NANOG; Nanog
Species	Human
GeneID (Human)	79923
GeneID	79923
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 153-305 of human Nanog (NP_079141.2)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9H9S0
Protein Weight	35kDa
Shipping	Ice bag