

Anti-Nanog Rabbit Polyclonal Antibody For Wb, Ihc-P, Elisa - Q9H9S0

Item Number: CM0026595



Introduction

Nanog Rabbit pAb is a polyclonal antibody validated for WB, IHC-P, and ELISA, with reactivity in human, mouse, and rat. Targets the NANOG transcription factor key to stem cell pluripotency and development.

[Learn More](#)

Parameter	Value
Product Name	Nanog Rabbit pAb
Remarks/Alias	NANOG; Nanog
Species	Human
GeneID (Human)	79923
GeneID	79923
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-105 of human Nanog (NP_079141.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9H9S0
Protein Weight	35kDa
Shipping	Ice bag