

Anti-Nanp Rabbit Polyclonal Antibody For Wb And Elisa - Q8Tbe9

Item Number: CM0017727



Introduction

Rabbit polyclonal antibody targeting human N-acylneuraminate-9-phosphatase (NANP/HDHD4). Verified for WB and ELISA, cross-reacts with mouse and rat. Key tool for sialic acid metabolism research. 28 kDa.

[Learn More](#)

Parameter	Value
Product Name	NANP Rabbit pAb
Remarks/Alias	HDHD4; C20orf147; dj694B14.3; NANP
Species	Human
GeneID (Human)	140838
GeneID	140838
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-248 of human NANP (NP_689880.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q8TBE9
Protein Weight	28kDa
Shipping	Ice bag