

Anti-Nans Rabbit Polyclonal Antibody For Wb, Elisa - Q9Nr45

Item Number: CM0028766



Introduction

Rabbit polyclonal antibody against human NANS, raised against a recombinant protein corresponding to amino acids 260-359. Suitable for Western blot and ELISA applications. Ideal for studying sialic acid metabolism.

[Learn More](#)

Parameter	Value
Product Name	NANS Rabbit pAb
Remarks/Alias	SAS; SEMDG; SEMDCG; HEL-S-100; NANS
Species	Human
GeneID (Human)	54187
GeneID	54187
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 260-359 of human NANS (NP_061819.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	Q9NR45
Protein Weight	40kDa
Shipping	Ice bag