

## Anti-Nans Rabbit Polyclonal Antibody For Wb, Elisa - Q9Nr45

Item Number: CM0021266



### Introduction

Rabbit polyclonal antibody against human NANS, validated for knockout studies. Suitable for western blotting and ELISA. NANS is a key enzyme in sialic acid biosynthesis, essential for brain and skeletal development.

[Learn More](#)

| Parameter        | Value                                                                                                             |
|------------------|-------------------------------------------------------------------------------------------------------------------|
| Product Name     | [KO Validated] NANS Rabbit pAb                                                                                    |
| Remarks/Alias    | SAS; SEMDG; SEMDCG; HEL-S-100; NS                                                                                 |
| Species          | Human                                                                                                             |
| GeneID (Human)   | 54187                                                                                                             |
| GeneID           | 54187                                                                                                             |
| Immunogen        | Recombinant fusion protein containing a sequence corresponding to amino acids 260-359 of human NANS (NP_061819.2) |
| Source           | Rabbit                                                                                                            |
| Category         | Polyclonal Antibodies                                                                                             |
| Application      | WB, ELISA                                                                                                         |
| Cross Reactivity | Human                                                                                                             |
| SWISS            | Q9NR45                                                                                                            |
| Protein Weight   | 40kDa                                                                                                             |
| Shipping         | Ice bag                                                                                                           |