

Anti-Nap1L5 Rabbit Polyclonal Antibody For Wb And Elisa - Q96Nt1

Item Number: CM0023690

Introduction

NAP1L5 Rabbit Polyclonal Antibody validated for WB and ELISA, reactive with human and mouse. Ideal for investigating histone chaperone function, nucleosome assembly, and brain-specific chromatin dynamics.

[Learn More](#)

Parameter	Value
Product Name	NAP1L5 Rabbit pAb
Remarks/Alias	DRLM; NAP1L5
Species	Human
GeneID (Human)	266812
GeneID	266812
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 30-120 of human NAP1L5 (NP_715638.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q96NT1
Protein Weight	20kDa
Shipping	Ice bag