

Anti-Napb Rabbit Polyclonal Antibody For Wb, Elisa - Q9H115

Item Number: CM0022761

Introduction

Rabbit polyclonal antibody against human NAPB (SNAP-BETA), validated for WB and ELISA. Suitable for studying ER-Golgi vesicular transport and related neurodevelopmental disorders in human, mouse, and rat samples.

[Learn More](#)

Parameter	Value
Product Name	NAPB Rabbit pAb
Remarks/Alias	SNAPB; DEE107; SNAP-BETA; NAPB
Species	Human
GeneID	63908
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-50 of human NAPB (NP_071363.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9H115
Protein Weight	34kDa
Shipping	Ice bag