

Anti-Napg (Gammاسnap) Rabbit Polyclonal Antibody For Wb And Elisa - Q99747

Item Number: CM0016854

Introduction

Rabbit polyclonal antibody against NAPG (GAMMASNAP), validated for WB and ELISA in human, mouse, and rat samples. Ideal for studying ER-to-Golgi vesicular transport.

[Learn More](#)

Parameter	Value
Product Name	NAPG Rabbit pAb
Remarks/Alias	GAMMASNAP; NAPG
Species	Human
Gene ID	8774
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-312 of human NAPG (NP_003817.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q99747
Protein Weight	35kDa
Shipping	Ice bag