

Anti-Napg Rabbit Polyclonal Antibody For Wb And Elisa - Q99747

Item Number: CM0027963

Introduction

NAPG Rabbit pAb is a polyclonal antibody against gamma-soluble NSF attachment protein, validated for WB and ELISA, with cross-reactivity in Human, Mouse, Rat. Suitable for intracellular trafficking research.

[Learn More](#)

Parameter	Value
Product Name	NAPG Rabbit pAb
Remarks/Alias	GAMMASNAP; NAPG
Species	Human
GeneID (Human)	8774
GeneID	8774
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 1-312 of human NAPG (NP_003817.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q99747
Protein Weight	35kDa
Shipping	Ice bag