

Anti-Naprt Polyclonal Antibody For Wb, Elisa - Q6Xqn6

Item Number: CM0018532



Introduction

Anti-NAPRT rabbit polyclonal antibody validated for WB and ELISA; detects human, mouse, rat NAPRT, a ~58 kDa cytosolic enzyme catalyzing the first step of NAD biosynthesis from nicotinic acid.

[Learn More](#)

Parameter	Value
Product Name	NAPRT Rabbit pAb
Remarks/Alias	NAPRT1; PP3856; NAPRT
Species	Human
Gene ID (Human)	93100
Gene ID	93100
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 309-538 of human NAPRT1 (NP_660202.3)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q6XQN6
Protein Weight	58kDa
Shipping	Ice bag