

Anti-Napsa Polyclonal Antibody For Wb, Elisa - 096009

Item Number: CM0019652

Introduction

NAPSA Rabbit pAb detects napsin-A, an aspartyl protease involved in surfactant processing. Suitable for WB and ELISA, with rat cross-reactivity. Ideal for lung and kidney research.

[Learn More](#)

Parameter	Value
Product Name	NAPSA Rabbit pAb
Remarks/Alias	NAPSA; KAP; Kdap; NAP1; NAPA; SNAPA; napsin-A
Species	Human
GeneID	9476
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 64-190 of human NAPSA (NP_004842.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Rat
SWISS	096009
Protein Weight	45kDa
Shipping	Ice bag