

Anti-Nars Rabbit Polyclonal Antibody For Wb, Elisa - O43776

Item Number: CM0013639



Introduction

NARS Rabbit Polyclonal Antibody validated for WB and ELISA against human and mouse asparagine--tRNA ligase. Targets a 63kDa cytoplasmic enzyme crucial for protein synthesis and cell signaling.

[Learn More](#)

Parameter	Value
Product Name	NARS Rabbit pAb
Remarks/Alias	NARS; ASNRS; NEDMILG; NEDMILEG
Species	Human
Gene ID (Human)	4677
Gene ID	4677
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 309-548 of human NARS (NP_004530.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	O43776
Protein Weight	63kDa
Shipping	Ice bag