

Anti-Nasp Polyclonal Antibody For Wb, Ihc-P, Elisa - P49321

Item Number: CM0027878



Introduction

Rabbit polyclonal antibody against human NASP (nuclear autoantigenic sperm protein) validated for WB, IHC-P, and ELISA; cross-reacts with mouse and rat. Suitable for studying histone chaperone networks and DNA replication.

[Learn More](#)

Parameter	Value
Product Name	NASP Rabbit pAb
Remarks/Alias	HMDRA1; FLB7527; PRO1999; NASP
Species	Human
Gene ID (Human)	4678
Gene ID	4678
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 180-449 of human NASP (NP_689511.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P49321
Protein Weight	85kDa
Shipping	Ice bag