


Anti-Nat1 Rabbit Monoclonal Antibody For If/Icc, Elisa - P18440

Item Number: CM0002913



Introduction

High-specificity anti-NAT1 rabbit monoclonal antibody validated for IF/ICC and ELISA. Cross-reacts with human, mouse, and rat NAT1 (34 kDa). Ideal for drug metabolism, pharmacogenomics, and cancer research studies.

[Learn More](#)

Parameter	Value
Product Name	NAT1 Rabbit mAb
Remarks/Alias	AAC1; MNAT; NATI; NAT-1; NAT1
Species	Human
GeneID (Human)	9
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 1-100 of human NAT1 (P18440)
Source	Rabbit
Category	Monoclonal Antibodies
Application	IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P18440
Protein Weight	34kDa
Shipping	Ice bag