

Anti-Nat1 Rabbit Polyclonal Antibody For Wb And Elisa - P18440

Item Number: CM0022140



Introduction

Rabbit polyclonal antibody targeting human NAT1 (arylamine N-acetyltransferase 1), validated for WB and ELISA. Recognizes a 34 kDa cytoplasmic enzyme involved in drug detoxification and carcinogen bioactivation, suitable for pharmacogenetic and toxicology research.

[Learn More](#)

Parameter	Value
Product Name	NAT1 Rabbit pAb
Remarks/Alias	AAC1; MNAT; NATI; NAT-1; NAT1
Species	Human
GeneID (Human)	9
GeneID	9
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-290 of human NAT1 (NP_001153646.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	P18440
Protein Weight	34kDa
Shipping	Ice bag