

Anti-Nat10 Rabbit Monoclonal Antibody () For Intracellular Flow Cytometry - Q9H0A0

Item Number: CM0011198



Introduction

rabbit monoclonal antibody against human NAT10, validated for intracellular flow cytometry. Targets the RNA acetyltransferase involved in mRNA stability and ribosome biogenesis. Ideal for IVD research and diagnostic development.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Human NAT10 mAb
Remarks/Alias	ALP; Kre33; NET43
Species	Human
GeneID (Human)	55226
GeneID	55226
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 750-1025 of human NAT10 (Q9H0A0)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC (intra)
Cross Reactivity	Human
SWISS	Q9H0A0
Protein Weight	116 kDa
Shipping	Ice bag