


Anti-Nat10 Rabbit Monoclonal Antibody For Wb, Ihc-P, Elisa - Q9H0A0

Item Number: CM0003145



Introduction

High-quality rabbit monoclonal antibody against NAT10 (RNA cytidine acetyltransferase). Validated for WB, IHC-P, ELISA; cross-reacts with human, mouse, rat. Ideal for studying mRNA acetylation and rRNA processing.

[Learn More](#)

Parameter	Value
Product Name	NAT10 Rabbit mAb
Remarks/Alias	ALP; Kre33; NET43; NAT10
Species	Human
GeneID (Human)	55226
GeneID	55226
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 750-1025 of human NAT10 (Q9H0A0)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9H0A0
Protein Weight	116kDa
Shipping	Ice bag