


Anti-Nat10 Rabbit Monoclonal Antibody [Kd Validated] For Wb, Ip, If/Icc, Ihc-P, Elisa - Q9H0A0

Item Number: CM0004693



Introduction

KD-validated rabbit monoclonal antibody targeting human NAT10 (RNA cytidine acetyltransferase). Suitable for WB, IP, IF/ICC, IHC-P, ELISA. Cross-reacts with mouse and rat. Ideal for studying RNA ac4C modifications and ribosome biogenesis.

[Learn More](#)

Parameter	Value
Product Name	Anti-NAT10 Rabbit Monoclonal Antibody [KD Validated]
Remarks/Alias	ALP; Kre33; NET43
Species	Human
GeneID (Human)	55226
GeneID	55226
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 676-1025 of human NAT10 (NP_078938.2).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IP, IF/ICC, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9H0A0
Protein Weight	116 kDa/107 kDa
Shipping	Ice bag