

Anti-Nat8 Rabbit Polyclonal Antibody For Wb, Ihc-P, Elisa - Q9Uhe5

Item Number: CM0030673

Introduction

High-quality rabbit polyclonal antibody targeting human NAT8, suitable for WB, IHC-P, and ELISA applications. Demonstrates cross-reactivity with mouse and rat. Ideal for studies on protein acetylation, apoptosis, and detoxification.

[Learn More](#)

Parameter	Specification
Product Name	NAT8 Rabbit pAb
Remarks/Alias	GLA, CML1, CCNAT, Hcml1, ATase2, TSC501, TSC510, NAT8
Species	Human
Gene ID	9027 (Human)
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 61-220 of human NAT8 (Q9UHE5)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9UHE5
Protein Weight	26kDa
Shipping	Ice bag