

Anti-Nat8B Rabbit Polyclonal Antibody (Cml2, Atase1) For Wb, Ihc-P, Elisa - Q9Uhf3

Item Number: CM0028096

Introduction

High-quality anti-NAT8B rabbit polyclonal antibody validated for WB, IHC-P, ELISA. Recognizes human, mouse, rat NAT8B (CML2/ATase1). Ideal for studying ER lysine acetylation and neurodegeneration research.

[Learn More](#)

Parameter	Value
Product Name	NAT8B Rabbit pAb
Remarks/Alias	CML2; Hcml2; NAT8BP; NAT8B
Species	Human
GeneID	51471
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 75-165 of human NAT8B (NP_057431.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9UHF3
Protein Weight	25kDa
Shipping	Ice bag