

Anti-Nat9 Rabbit Polyclonal Antibody For Wb And Elisa - Q9Bte0

Item Number: CM0030690

Introduction

High-quality rabbit polyclonal antibody against human NAT9 (alpha/beta-tubulin-N-acetyltransferase 9), suitable for Western blot and ELISA. Ideal for tubulin acetylation and neurobiology research.

[Learn More](#)

Parameter	Value
Product Name	NAT9 Rabbit pAb
Remarks/Alias	EBSP; hNATL; NAT9
Species	Human
GenelD (Human)	26151
GenelD	26151
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-207 of human NAT9 (NP_056469.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	Q9BTE0
Protein Weight	23kDa
Shipping	Ice bag