


Anti-Natriuretic Peptides A Rabbit Monoclonal Antibody For Wb, If-P, Ihc-P, Elisa - P01160

Item Number: CM0007839



Introduction

Anti-Natriuretic peptides A rabbit monoclonal antibody validated for WB, IF-P, IHC-P, ELISA; cross-reacts with human, mouse, rat.

Recombinant immunogen targets the cardiac hormone regulating blood pressure, fluid balance, and cardiovascular remodeling. P01160.

[Learn More](#)

Parameter	Value
Product Name	Natriuretic peptides A Rabbit mAb
Remarks/Alias	ANF; ANP; CDD; CDP; PND; ATFB6; ATRST2; CDD-ANF
Species	Human
GeneID (Human)	4878
GeneID	4878
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 26-151 of human Natriuretic peptides A (NP_006163.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF-P, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P01160
Protein Weight	16 kDa
Shipping	Ice bag