


Anti-Nav1.5/Scn5A Monoclonal Antibody For Wb, Ihc-P, Elisa - Q14524

Item Number: CM0009314



Introduction

High-quality Anti-Nav1.5/SCN5A Rabbit Monoclonal Antibody validated for WB, IHC-P, ELISA. Specific for human Nav1.5, crucial for cardiac electrophysiology research and IVD assay development.

[Learn More](#)

Parameter	Value
Product Name	Nav1.5/SCN5A Rabbit mAb
Remarks/Alias	HB1; HB2; HH1; IVF; VF1; HBBD; ICCD; LQT3; SSS1; CDCD2; CMD1E; CMPD2; PFHB1; Nav1.5
Species	Human
Gene ID (Human)	6331
Gene ID	6331
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 273-357 of human Nav1.5/SCN5A (NP_000326.2)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human
SWISS	Q14524
Protein Weight	227 kDa/223 kDa/221 kDa/225 kDa
Shipping	Ice bag