

Anti-Naxe Polyclonal Antibody For Wb, Elisa - Q8Ncw5

Item Number: CM0015273

Introduction

Rabbit polyclonal antibody targeting human NAXE (AIBP/APOA1BP). Detects human, mouse, and rat NAXE. Suitable for WB and ELISA. Targets a mitochondrial NAD(P)H-hydrate epimerase involved in cholesterol efflux and vascular regulation.

[Learn More](#)

Parameter	Value
Product Name	NAXE Rabbit pAb
Remarks/Alias	AIBP; PEBEL; YJEFN1; APOA1BP; NAXE
Species	Human
GeneID (Human)	128240
GeneID	128240
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 48-170 of human NAXE (NP_658985.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q8NCW5
Protein Weight	32kDa
Shipping	Ice bag