

Anti-Nbas Polyclonal Antibody For Wb, If/Icc, Elisa - A2Rrp1

Item Number: CM0025204

Introduction

Affinity-purified rabbit anti-NBAS polyclonal antibody validated for Western blot, IF/ICC, and ELISA. Reacts with human and rat NBAS. A research tool for Golgi-to-ER retrograde transport and neuroblastoma studies.

[Learn More](#)

Parameter	Value
Product Name	NBAS Rabbit pAb
Remarks/Alias	NAG; SOPH; ILFS2; NBAS
Species	Human
Gene ID (Human)	51594
Gene ID	51594
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 1-70 of human NBAS (NP_056993.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Rat
SWISS	A2RRP1
Protein Weight	269kDa
Shipping	Ice bag