


Anti-Nbr1 Rabbit Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - Q14596

Item Number: CM0030767



Introduction

Anti-NBR1 rabbit polyclonal antibody validated for WB, IHC-P, IF/ICC, and ELISA with human, mouse, rat reactivity. Autophagy adapter involved in innate immunity and protein aggregation. Ideal for autophagy and cancer research.

[Learn More](#)

Parameter	Value
Product Name	NBR1 Rabbit pAb
Remarks/Alias	IAI3B; M17S2; MIG19; 1A1-3B; NBR1
Species	Human
GeneID (Human)	4077
GeneID	4077
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-85 of human NBR1 (NP_005890.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q14596
Protein Weight	107kDa
Shipping	Ice bag