

Anti-Nbs1/Nbn Monoclonal Antibody For Wb, Ip, Elisa - O60934

Item Number: CM0004434



Introduction

Rabbit monoclonal antibody against human NBS1/NBN, validated for WB, IP, ELISA. Nibrin is a component of the MRN complex essential for DNA double-strand break repair and telomere maintenance.

[Learn More](#)

Parameter	Value
Product Name	NBS1/NBN Rabbit mAb
Remarks/Alias	ATV; NBS; P95; NBS1; AT-V1; AT-V2; NBS1/NBN
Species	Human
GeneID (Human)	4683
GeneID	4683
Immunogen	Synthetic Peptide A synthetic peptide corresponding to a sequence within amino acids 655-754 of human NBS1/NBN (O60934).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IP, ELISA
Cross Reactivity	Human
SWISS	O60934
Protein Weight	85kDa
Shipping	Ice bag