

Anti-Nbs1/Nbn Polyclonal Antibody For Wb, If/Icc, Elisa - O60934

Item Number: CM0027592



Introduction

CamelBio offers a rabbit polyclonal anti-NBS1/NBN (Nibrin) antibody, validated for WB, IF/ICC, and ELISA, with cross-reactivity to human, mouse, and rat. Ideal for DNA repair and telomere integrity studies.

[Learn More](#)

Parameter	Value
Product Name	NBS1/NBN Rabbit pAb
Remarks/Alias	ATV, NBS, P95, NBS1, AT-V1, AT-V2, NBS1/NBN
Species	Human
Gene ID	4683
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 650-750 of human NBS1/NBN (NP_002476.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	O60934
Protein Weight	85kDa
Shipping	Ice bag