

Anti-Nbs1/Nbn Rabbit Polyclonal Antibody (Ko Validated) For Wb, Ihc-P, Elisa - O60934

Item Number: CM0017899



Introduction

KO-validated rabbit polyclonal antibody against NBS1/NBN (Nibrin). Suitable for WB, IHC-P, ELISA. Recognizes human, mouse, rat. Targets DNA repair, telomere maintenance. SWISS: O60934.

[Learn More](#)

Parameter	Value
Product Name	[KO Validated] NBS1/NBN Rabbit pAb
Remarks/Alias	ATV, NBS, P95, NBS1, AT-V1, AT-V2, BN
Species	Human
GeneID (Human)	4683
GeneID	4683
Immunogen	Recombinant protein / Recombinant fusion protein containing a sequence corresponding to amino acids 520-754 of human NBS1/NBN (NP_002476.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	O60934
Protein Weight	85kDa
Shipping	Ice bag