

Anti-Ncald Rabbit Polyclonal Antibody For Wb And Elisa - P61601

Item Number: CM0031026

Introduction

Rabbit polyclonal antibody targeting human neurocalcin-delta (NCALD). Suitable for WB and ELISA applications, with cross-reactivity to mouse and rat. Ideal for calcium signaling and neuronal tissue research.

[Learn More](#)

Parameter	Value
Product Name	NCALD Rabbit pAb
Remarks/Alias	NCALD
Species	Human
GeneID (Human)	83988
GeneID	83988
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-193 of human NCALD (NP_114430.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P61601
Protein Weight	22kDa
Shipping	Ice bag