

Anti-Ncam Rabbit Monoclonal Antibody For Wb, Ihc-P, Elisa - P13591

Item Number: CM0003959



Introduction

Rabbit monoclonal antibody targeting human NCAM (CD56, NCAM1), validated for Western blot, immunohistochemistry, and ELISA. Ideal for studying neuron adhesion, neurite outgrowth, and virus entry mechanisms.

[Learn More](#)

Parameter	Value
Product Name	NCAM Rabbit mAb
Remarks/Alias	CD56; NCAM; MSK39
Species	Human
GeneID (Human)	4684
GeneID	4684
Immunogen	Synthetic Peptide A synthetic peptide corresponding to a sequence within amino acids 758-858 of human NCAM (NP_000606.3).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human
SWISS	P13591
Protein Weight	95kDa
Shipping	Ice bag