

Anti-Ncam1/Cd56 Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - P13591

Item Number: CM0030658



Introduction

Rabbit polyclonal antibody against human NCAM1/CD56, validated for WB, IHC-P, IF/ICC, and ELISA. Cross-reacts with mouse and rat. Essential for neural adhesion and neurobiology research.

[Learn More](#)

Parameter	Value
Product Name	NCAM1 / CD56 Rabbit pAb
Remarks/Alias	CD56; NCAM; MSK39; NCAM1 / CD56
Species	Human
GeneID (Human)	4684
GeneID	4684
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 759-858 of human NCAM1 / CD56 (NP_851996.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P13591
Protein Weight	95kDa
Shipping	Ice bag