

Anti-Ncam1/Cd56 Rabbit Polyclonal Antibody For Wb, Ihc-P, Elisa - P13591

Item Number: CM0013084



Introduction

High-quality rabbit polyclonal antibody against human NCAM1/CD56, validated for WB, IHC-P, and ELISA. Recognizes mouse and rat homologs; suitable for studies of neural cell adhesion, neurite outgrowth, and virus receptor interactions.

[Learn More](#)

Parameter	Specification
Product Name	NCAM1 / CD56 Rabbit pAb
Remarks/Alias	CD56; NCAM; MSK39; NCAM1 / CD56
Species	Human
Gene ID (Human)	4684
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 759-858 of human NCAM1 / CD56 (NP_851996.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Mouse, Rat
SWISS	P13591
Protein Weight	95kDa
Shipping	Ice bag