

## Anti-Ncam/Cd56 Rabbit Monoclonal Antibody For Wb, If-P, Ihc-P, Elisa - P13595

Item Number: CM0005964



### Introduction

CamelBio's NCAM/CD56 Rabbit Monoclonal Antibody (P13595) validated for WB, IF-P, IHC-P, ELISA. Recognizes mouse and rat NCAM isoforms (67-119 kDa). Ideal for neuroscience and cell adhesion studies.

[Learn More](#)

| Parameter        | Value                                                                                                                    |
|------------------|--------------------------------------------------------------------------------------------------------------------------|
| Product Name     | NCAM /CD56 Rabbit mAb                                                                                                    |
| Remarks/Alias    | CD56; Ncam; E-NCAM; NCAM-1; NCAM /CD56                                                                                   |
| Species          | Mouse                                                                                                                    |
| GeneID (Human)   | 17967                                                                                                                    |
| GeneID           | 17967                                                                                                                    |
| Immunogen        | Recombinant fusion protein containing a sequence corresponding to amino acids 20-711 of mouse NCAM /CD56(NP_001106675.1) |
| Source           | Rabbit                                                                                                                   |
| Category         | Monoclonal Antibodies                                                                                                    |
| Application      | WB, IF-P, IHC-P, ELISA                                                                                                   |
| Cross Reactivity | Mouse, Rat                                                                                                               |
| SWISS            | P13595                                                                                                                   |
| Protein Weight   | 67kDa/80kDa/93kDa/119kDa                                                                                                 |
| Shipping         | Ice bag                                                                                                                  |