

Anti-Ncaph Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - Q15003

Item Number: CM0025483



Introduction

Rabbit polyclonal antibody targeting NCAPH (condensin complex subunit 2) for WB, IF/ICC, ELISA applications. Cross-reacts with human. Ideal for chromosome condensation and cell cycle research. SWISS: Q15003.

[Learn More](#)

Parameter	Value
Product Name	NCAPH Rabbit pAb
Remarks/Alias	CAPH; BRRN1; CAP-H; MCPH23; NCAPH
Species	Human
Gene ID (Human)	23397
Gene ID	23397
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-250 of human NCAPH (NP_056156.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human
SWISS	Q15003
Protein Weight	83kDa
Shipping	Ice bag