


Anti-Ncbp2 Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - P52298

Item Number: CM0027239



Introduction

NCBP2 Rabbit Polyclonal Antibody validated for Western blot, immunofluorescence/ICC, and ELISA. Detects human and mouse NCBP2 (CBP20), a key subunit of the nuclear cap-binding complex essential for pre-mRNA splicing, export, and translation regulation.

[Learn More](#)

Parameter	Value
Product Name	NCBP2 Rabbit pAb
Remarks/Alias	CBC2; NIP1; CBP20; PIG55; NCBP2
Species	Human
GeneID (Human)	22916
GeneID	22916
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-156 of human NCBP2 (NP_031388.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse
SWISS	P52298
Protein Weight	18kDa
Shipping	Ice bag