

Anti-Ncf4/P40-Phox Monoclonal Antibody For Wb, Ihc-P, Ip, Elisa - Q15080

Item Number: CM0003235



Introduction

Rabbit monoclonal antibody against human NCF4/p40-phox (Q15080), validated for WB, IHC-P, IP, ELISA. Suitable for studying phagocyte NADPH oxidase function and chronic granulomatous disease.

[Learn More](#)

Parameter	Value
Product Name	NCF4/p40-phox Rabbit mAb
Remarks/Alias	NCF; CGD3; P40PHOX; SH3PXD4; NCF4/p40-phox
Species	Human
GeneID	4689
Immunogen	Synthetic peptide corresponding to amino acids 240-339 of human NCF4/p40-phox (Q15080)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, IP, ELISA
Cross Reactivity	Human
SWISS	Q15080
Protein Weight	39kDa
Shipping	Ice bag