

Anti-Ncf4/P40-Phox Polyclonal Antibody For Wb, Elisa - Q15080

Item Number: CM0024308



Introduction

Rabbit polyclonal antibody against human NCF4/p40-phox, validated for Western blot (WB) and ELISA. Detects endogenous p40-phox in human, mouse, and rat samples. Ideal for studying NADPH oxidase complex in phagocytic cells.

[Learn More](#)

Parameter	Value
Product Name	NCF4/p40-phox Rabbit pAb
Remarks/Alias	NCF; CGD3; P40PHOX; SH3PXD4; NCF4/p40-phox
Species	Human
GeneID	4689
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 1-190 of human NCF4/p40-phox (NP_038202.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q15080
Protein Weight	39kDa
Shipping	Ice bag