

Anti-Nck1 Monoclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - P16333

Item Number: CM0010320



Introduction

NCK1 rabbit monoclonal antibody validated for WB, IHC-P, IF/ICC, ELISA applications. Detects human, mouse, rat NCK1 (43kDa, P16333). Targets SH2/SH3 adapter protein involved in growth factor signaling, translation regulation, and DNA damage response.

[Learn More](#)

Parameter	Value
Product Name	NCK1 Rabbit mAb
Remarks/Alias	NCK; nck-1; NCKalpha; NCK1
Species	Human
GeneID (Human)	4690
GeneID	4690
Immunogen	Synthetic Peptide (A synthetic peptide corresponding to a sequence within amino acids 200-300 of human NCK1 (P16333))
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P16333
Protein Weight	43kDa
Shipping	Ice bag