

Anti-Nck1 Rabbit Polyclonal Antibody For Wb And Elisa - P16333

Item Number: CM0018253



Introduction

Rabbit polyclonal antibody targeting NCK1 (NCKalpha), an SH2/SH3 adapter protein. Suitable for WB and ELISA; cross-reacts with human and mouse. For investigating receptor signaling pathways and DNA damage response.

[Learn More](#)

Parameter	Value
Product Name	NCK1 Rabbit pAb
Remarks/Alias	NCK; nck-1; NCKalpha; NCK1
Species	Human
Gene ID (Human)	4690
Gene ID	4690
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 58-377 of human NCK1 (NP_006144.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	P16333
Protein Weight	43kDa
Shipping	Ice bag